



# Recombinant Rabbit C-reactive protein (CRP)

<b>Product Code</b>	CSB-EP005991RB-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P02742
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AGMHKKAFFV PKESDNSYVS LNAQLKKPLK AFTVCLYFYT DLSMTRGYSI FSYATRRQFN EILLFWSKDI GYSFVGGDE IIFKVSDIPV DPTHLCASWE SSTGIAELWV DGKPMVRKSL KKG YILGPEA SIILGQDQDS FGGSFEKQQS LVGDIGNVNM WDYALSPEEI NTIYAGGTFS PNVLDWRELT YQVRGEVHVK PQLWP
<b>Source</b>	E.coli
<b>Target Names</b>	CRP
<b>Protein Names</b>	Recommended name: C-reactive protein
<b>Expression Region</b>	21-225
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein belongs to the pentaxin family. It is involved in several host defense related functions based on its ability to recognize foreign pathogens and damaged cells of the host and to initiate their elimination by interacting with humoral and cellular effector systems in the blood. Consequently, the level of this protein in plasma increases greatly during acute phase response to tissue injury, infection, or other inflammatory stimuli.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.