



Recombinant Drosophila simulans Alpha-methyldopa hypersensitive protein (amd)

Product Code	CSB-MP305546DMJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81893
Product Type	Recombinant Protein
Immunogen Species	Drosophila simulans (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	TSTSYPSIVG EMLASGFGVI GFSWICSPAC TELEVVMDW LAKFLKPPAH FQHASDGGGG GVIQGSASEA VLVAVLAARE QAVASYRESH PELSESEVRG RLVAYSSDQS NSCIEKAGVL AAMPIRLLPA GEDFVLRGDT LRGAIEEDVA AGRIPVICVA TLGTTGTCAV DDIESLSAVC EEFKVWLHVD AAYAGGAFAL EECSDLRKGL DRVDSLNFNL HKFMLVNFDC SAMWLRDANK VVDSFNVDRV YLKHKHEGQS QIPDFRHWQI PLGRRFRALK VWITFRTLGA EGLRNHVRKH IELAKQFEQL VLKDSRFELV APSALGLV
Source	Mammalian cell
Target Names	amd
Protein Names	Recommended name: Alpha-methyldopa hypersensitive protein EC= 4.1.1.-
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.