



Recombinant *Drosophila melanogaster* cAMP-dependent protein kinase type II regulatory subunit (Pka-R2)

Product Code	CSB-MP305629DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81900
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSSDSSRRIQ VPEELKEVLL QFSISFLVEQ PPDVIDYAVE YFTKLQSERP SVSHTDQSTD DQLSVNSQDA DAEPVMASS RRKSVFAEAY DPEADDDDDG ATAVFPKTDE QRARLVESVK NVLLFRSLEK EQMNQVLDAM FERKVQPGDF IIRQGDDGDN FYVIESGVYK VYINDKHINT YNHTGLFGEL ALLYNMPRAA TVQAETSGLL WAMDRQTFRR ILLKSAFRKR KMYEELLNSV PMLKALQNYE RMNLADALVS KSYDNGERII KQGDAADGMY FIEEGTVSVR MDQDDAEVEI SQLGKGQYFG ELALVTHRPR AASVYATGGV VKLAFLDVKA FERLLGPCMD IMKRNIDDYE SQLVKIFGSK NNITDTR
Source	Mammalian cell
Target Names	Pka-R2
Protein Names	Recommended name: cAMP-dependent protein kinase type II regulatory subunit
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.