



Recombinant Arabidopsis thaliana Thylakoid lumenal 17.4 kDa protein, chloroplastic (At5g53490)

Product Code	CSB-EP305777DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81760
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	ANQ RLPPLSTEPD RCEKAFVGNT IGQANGVYDK PLDLRFCDYT NDQTNLKGKT LSAALMVGAK FDGADMTEVV MSKAYAVEAS FKGVNFTNAV IDRVNFGKSN LKGAVFRNTV LSGSTFEEAN LEDVVFEDTI IGYIDLQKIC RNESINEEGR LVLGCR
Source	E.coli
Target Names	TL17
Protein Names	Recommended name: Thylakoid lumenal 17.4 kDa protein, chloroplastic Alternative name(s): P17.4
Expression Region	78-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.