



Recombinant *Geobacillus stearothermophilus* 30S ribosomal protein S2 (rpsB)

Product Code	CSB-BP305831GFM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81289
Product Type	Recombinant Protein
Immunogen Species	<i>Geobacillus stearothermophilus</i> (<i>Bacillus stearothermophilus</i>)
Purity	>85% (SDS-PAGE)
Sequence	SVISMKQLL EAGVHFGWQT RRWNPKMKKY IFTERNGIYL IDIQKTVKKV EEAYNFVREL AANGGKILFV GTKRQAQESV KEEAERCGMF YVNQRRIGGT LTNFATLQKR IKRLREIEKM EEDGVFDVLP KKEVIGLKKE KERLEKFIGG IKDMKELPDA LFVIDPRKER IAVAEARKLN IPIIGIVDTN NDPCEGGYVI PANDDAIRAV KLLTSKIADA VLEAKQGEEA AVAAE
Source	Baculovirus
Target Names	rpsB
Protein Names	Recommended name: 30S ribosomal protein S2 Alternative name(s): BS2a
Expression Region	2-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.