




# Recombinant Human Cellular tumor antigen p53 (TP53)

<b>Product Code</b>	CSB-YP024077HU
<b>Relevance</b>	Acts as a tumor suppressor in many tumor types; induces growth arrest or apoptosis depending on the physiological circumstances and cell type. Involved in cell cycle regulation as a trans-activator that acts to negatively regulate cell division by controlling a set of genes required for this process. One of the activated genes is an inhibitor of cyclin-dependent kinases. Apoptosis induction seems to be mediated either by stimulation of BAX and FAS antigen expression, or by repression of Bcl-2 expression. In cooperation with mitochondrial PPIF is involved in activating oxidative stress-induced necrosis; the function is largely independent of transcription. Induces the transcription of long intergenic non-coding RNA p21 (lincRNA-p21) and lincRNA-Mkl1. LincRNA-p21 participates in TP53-dependent transcriptional repression leading to apoptosis and seem to have to effect on cell-cycle regulation. Implicated in Notch signaling cross-over. Prevents CDK7 kinase activity when associated to CAK complex in response to DNA damage, thus stopping cell cycle progression. Isoform 2 enhances the transactivation activity of isoform 1 from some but not all TP53-inducible promoters. Isoform 4 suppresses transactivation activity and impairs growth suppression mediated by isoform 1. Isoform 7 inhibits isoform 1-mediated apoptosis. Regulates the circadian clock by repressing CLOCK-ARNTL/BMAL1-mediated transcriptional activation of PER2
<b>Abbreviation</b>	Recombinant Human TP53 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P04637
<b>Alias</b>	Antigen NY-CO-13 Phosphoprotein p53 Tumor suppressor p53
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MEEPQSDPSVEPPLSQETFSDLWKLLPENNVLSPLPSQAMDDLMLSPDDIEQ WFTEDPGPDEAPRMPEAAPPVAPAPAAPTTPAAPAPAPSWPLSSSVPSQKTYQ GSYGFRGLGFLHSGTAKSVTCTYSPALNKMFCQLAKTCPVQLWVDSTPPPGTR VRAMAIYKQSQHMTEVVRRCPPHHERCSDSDGLAPPQHILIRVEGNLRVEYLDD RNTFRHSVVVPYEPPEVGS DCTTIHYNMCMNSSCMGGMNRRPILTIITLEDSS GNLLGRNSFEVRVCACPGDRDRRTEENLRKKGEPHHELPPGSTKRALPNNTS SSPQPKKKPLDGEYFTLQIRGRERFEMFRELNEALELKDAQAGKEPGGSRAH SSHLKSKKGQSTSRHKKLMFKTEGPDSD
<b>Research Area</b>	Cell Biology



<b>Source</b>	Yeast
<b>Target Names</b>	TP53
<b>Protein Names</b>	Recommended name: Cellular tumor antigen p53 Alternative name(s): Antigen NY-CO-13 Phosphoprotein p53 Tumor suppressor p53
<b>Expression Region</b>	1-393aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	45.7kDa
<b>Protein Length</b>	Full Length

**Image**


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution** We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life** The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.