



Recombinant Rat Alpha-enolase (Eno1)

Product Code	CSB-YP007670RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04764
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	SILKIHARE IFDSRGNPTV EVDLYTAKGL FRAAVPSGAS TGIYEALELR DNDKTRFMGK GVS KAVEHIN KTIAPALVSK KLN VVEQEKI DQLMIEMDGT ENKSKFGANA ILGVSLAVCK AGAVEKGVPL YRHIADLAGN PEVILPVPF NVINGGSHAG NKLAMQEFMI LPVGASSFRE AMRIGAEVYH NLKNVIKEY GKDATNVGDE GGFAPNILEN KEALELLKSA IAKAGYTDQV VIGMDVAASE FYRAGKYDLD FKSPDDASRY ITPDQLADLY KSF IKDYPVV SIEDPFDQDD WDAWQKFTAT AGIQVVGDDL TVTNP KRIAK AAGEKSCNCL LLKVNQIGSV TESLQACKLA QSN GWGMVS HRSGETEDTF IADLVVGLCT GQIKTGAPCR SERLAKYNQI LRIEEEELGSK AKFAGRSFRN PLAK
Source	Yeast
Target Names	Eno1
Protein Names	Recommended name: Alpha-enolase EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase Enolase 1 Non-neural enolase Short name= NNE
Expression Region	2-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes one of three enolase isoenzymes found in mammals; it encodes alpha-enolase, a homodimeric soluble enzyme, and also encodes a shorter monomeric structural lens protein, tau-crystallin. The two proteins are made from the same message. The full length protein, the isoenzyme, is found in the cytoplasm. The shorter protein is produced from an alternative translation start, is localized to the nucleus, and has been found to bind to an element in the c-myc promoter. A pseudogene has been identified that is located on the other arm of the same chromosome.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.