



Recombinant Bovine Clathrin light chain A (CLTA)

Product Code	CSB-MP005591BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04973
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MAELDPFGVP AGGPALGNGV AGEEDPAAAF LAQQESEIAG IENDEAFAIL DGGAPGSQPH GEPPGIPDAV DGVTNGDYQ ESNGPTDSYA AISQVDRQLQS EPESIRKWRE EQTERLEALD ANSRKQEA EW KEKAIKELDE WYARQDEQLQ KTKANNRVAD EAFYKQPFAD VIGYVTNINH PCYSLEQAAE EAFVNDIEES SPGTEWERVA RLCDFNPKSS KQAKDVSRMR SVLISLQKAP LVH
Source	Mammalian cell
Target Names	CLTA
Protein Names	Recommended name: Clathrin light chain A Short name= Lca
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	Clathrin is a large, soluble protein composed of heavy and light chains. It functions as the main structural component of the lattice-type cytoplasmic face of coated pits and vesicles which entrap specific macromolecules during receptor-mediated endocytosis. This gene encodes one of two clathrin light chain proteins which are believed to function as regulatory elements. Alternative splicing results in multiple transcript variants.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.