



Recombinant Drosophila melanogaster Protein lethal (2)essential for life

Product Code	CSB-BP306768DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P82147
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSVVPLMFRD WWDELDFPMR TSRLLDQHFG QGLKRDDLMS SVWNSRPTVL RSGYLRPWHT NSLQKQESGS TLNIDSEKFE VILDVQQFSP SEITVKVADK FVIVEGKHEE KQDEHGYVSR QFSRRYQLPS DVNPDTVTSS LSSDGLLTIK APMKALPPPQ TERLVQITQT GPSSKEDNAK KVETSTA
Source	Baculovirus
Target Names	l(2)efl
Protein Names	Recommended name: Protein lethal(2)essential for life Alternative name(s): Protein Efl21
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.