



Recombinant Human Neuroendocrine protein 7B2 (SCG5), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP020810HU |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P05408 |
| Storage Buffer | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| Product Type | Recombinant Proteins |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | YSPRTPDRVSEADIQRLLHGVMEQLGIARPRVEYPAHQAMNLVGPQSIEGGA HEGLQHLGPFGNIPNIVAEELTGDNIPKDFSEDQGYDPPNPCPVGKTADDGCL ENTPDTAEFSSREFQLHQHLFDPEHDYPGLGKWNKKLLYEKMKGGE |
| Source | E.coli |
| Target Names | SCG5 |
| Protein Names | Recommended name: Neuroendocrine protein 7B2 Alternative name(s): Pituitary polypeptide Secretogranin V Secretogranin-5 Secretory granule endocrine protein I Cleaved into the following 2 chains: 1. N-terminal peptide 2. C-termi |
| Expression Region | 27-176aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |