



# Recombinant Human Serine protease hepsin (HPN), partial

<b>Product Code</b>	CSB-BP010704HU1
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P05981
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	RSDQEP LYPVQVSSAD ARLMVFDKTE GTWRLLCSSR SNARVAGLSC EEMGFLRALT HSELDVRTAG ANGTSGFFCV DEGRLPHTQR LLEVISVDCD PRGRFLAAIC QDCGRRKLPV DRIVGGRDTS LGRWPWQVSL RYDGAHLCCG SLLSGDWVLT AAHCFPERNR VLSRWRVFAG AVAQASPHGL QLGVQAVVYH GGYLPFRDPN SEENSDIAL VHLSSPLPLT EYIQPVCLPA AGQALVDGKI CTVTGWGNTQ YYGQQAGVLQ EARVPIISND VCNGADFYGN QIKPKMFCA G YPEGGIDACQ GDSGGPFVCE DSISRTPRWR LCGIVSWG TG CALAQKPGVY TKVSDFREWI FQAIKTHSEA SGMVTQL
<b>Source</b>	Baculovirus
<b>Target Names</b>	HPN
<b>Protein Names</b>	Recommended name: Serine protease hepsin EC= 3.4.21.106 Alternative name(s): Transmembrane protease serine 1 Cleaved into the following 2 chains: 1. Serine protease hepsin non-catalytic chain 2. Serine protease hepsin catalytic chain
<b>Expression Region</b>	45-417
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.