



Recombinant Arabidopsis thaliana Putative defensin-like protein 239 (SCRL17)

Product Code	CSB-MP307005DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P82636
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	GLIRECD YTMEGYGPCG PNGRSLCLDT FWKTPPSPGL SKNLEGCRCE DRGKRPIIFT GLSHTCRCCW SYGGGSDD
Source	Mammalian cell
Target Names	SCRL17
Protein Names	Recommended name: Putative defensin-like protein 239 Alternative name(s): Putative S locus cysteine-rich-like protein 17 Short name= Protein SCRL17 Short name= SCR-like protein 17
Expression Region	24-98
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.