



Recombinant Agkistrodon contortrix contortrix Thrombin-like enzyme contortrixobin

Product Code	CSB-MP307045AET
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P82981
Product Type	Recombinant Protein
Immunogen Species	Agkistrodon contortrix contortrix (Southern copperhead)
Purity	>85% (SDS-PAGE)
Sequence	VVGGDECNIN EHRFLVAIFN SNGFVCSGTL INQEWVLTAA HCDSTDFQIK LGAHKKVNLN EDEQIRNPKE KFICPNKKND EVLDKDIMLI KLDSRVSNSE HIAPLSLPSS PPSVGSVCHI MGWGSITPIE VTFPDVPHCA YINLLDDAAC QPGYPEVLPE YRTLCA GILE GGDTCNYDS GGPLICNGQF QGIVSYGAHP CGQSLKPGIY TKVFDYNDWI QSIAGNTAA TCPP
Source	Mammalian cell
Protein Names	Recommended name: Thrombin-like enzyme contortrixobin Short name= SVTLE xx EC= 3.4.21.- Alternative name(s): Fibrinogen-clotting enzyme Snake venom serine protease Short name= SVSP Venombin B
Expression Region	1-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.