



Recombinant Mouse Cathepsin L1 (CtSl1), partial

| | |
|--------------------------|---|
| Product Code | CSB-BP006193MO |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P06797 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | IPKSVDW REKGCVTPVK NQGQCGSCWA FSASGCLEGQ MFLKTGKLIS LSEQNLVDCS HAQGNQGCNG GLMDFAFQYI KENGGLDSEE SYPYEAKDGS CKYRAEFAVA NDTGFVDIPQ QEKALMKAVA TVGPISVAMD ASHPSLQFYS SGIYYEPNCS SKNLDHGVLL VGYGYEGT |
| Source | Baculovirus |
| Target Names | CtSl |
| Protein Names | Recommended name: Cathepsin L1 EC= 3.4.22.15 Alternative name(s): Major excreted protein Short name= MEP p39 cysteine proteinase Cleaved into the following 2 chains: 1. Cathepsin L1 heavy chain 2. Cathepsin L1 light chain |
| Expression Region | 114-288 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |