



# Recombinant Mouse NADP-dependent malic enzyme (Me1)

<b>Product Code</b>	CSB-BP013632MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P06801
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEPRAPRRRH THQRGYLLTR DPHLNKDLAF TLEERQQLNI HGLLPPCIIS QELQVLRRIK NFERLNSDFD RYLLMLDLQD RNEKLFYSVL MSDVEKFMPI VYTPTVGLAC QQYSLAFRKP RGLFISIHDK GHIASVLNAW PEDVVKAIVV TDGERILGLG DLGCNGMGIP VGKLALYTAC GGVNPQQCLP ITLDVGTENE ELLKDPLYIG LRHRRVRGPE YDAFLDEFME AASSKYGMNC LIQFEDFANR NAFRLLNKYR NKYCTFNDDI QGTASVAVAG LLAALRITKN KLSDQTVLFQ GAGEAALGIA HLVVMAMEKE GLSKENARKK IWLVDKGLI VKGRASLTEE KEVFAHEHEE MKNLEAIVQK IKPTALIGVA AIGGAFTEQI LKDMAAFNER PIIFALSNT SKAECSAEQC YKVTKGRAIF ASGSPFDPVT LPDGRTLFPG QGNNSYVFPG VALGVVACGL RHIDDKVFLT TAEVISQQVS DKHLQEGRLY PPLNTIRGVS LKIAVKIVQD AYKEKMATVY PEPQNKEEFV SSQMYSTNYD QILPDCYPWP AEVQKIQTKV NQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	Me1
<b>Protein Names</b>	Recommended name: NADP-dependent malic enzyme Short name= NADP-ME EC= 1.1.1.40 Alternative name(s): Malic enzyme 1
<b>Expression Region</b>	1-572
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a cytosolic, NADP-dependent enzyme that generates NADPH for fatty acid biosynthesis. The activity of this enzyme, the reversible oxidative decarboxylation of malate, links the glycolytic and citric acid cycles. The regulation of expression for this gene is complex. Increased expression can result from elevated levels of thyroid hormones or by higher proportions of carbohydrates in the diet.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.