



Recombinant Chicken Calretinin (CALB2)

Product Code	CSB-EP004433CH-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P07090
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAGQRAPHLH LAELSASQFL DVWRHFDADG NGYIEGKELE NFFQELESAR KGTGVDSKRD SLGDKMKEFM HKYDKNADGK IEMAELAQIL PTEENFLLCF RQHVGSSSEF MEAWRRYDTD RSGYIEANEL KGFLSDLLKK ANRPYDEAKL QEYTQTILRM FDMNGDGKLG LSEMSRLLPV QENFLLKFQG MKLSSEEFNA IFAFYDKDGS GFIDEHELDA LLKDLYEKNK KEMSIQQLTN YRRSIMNLSLSD GGKLYRKELE VVLCSEPPL
Source	E.coli
Target Names	CALB2
Protein Names	Recommended name: Calretinin Short name= CR
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes an intracellular calcium-binding protein belonging to the troponin C superfamily. Members of this protein family have six EF-hand domains which bind calcium. Three alternatively spliced transcript variants that encode different proteins have been described.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.