



Recombinant Human Tyrosine 3-monoxygenase (TH)

Product Code	CSB-MP023470HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P07101
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	PTPDATTPQ AKGFRRVAVSE LDAQQAEAIM VRGQGAPGPS LTGSPWPGTA APAASYTTP RSPRFIGRRQ SLIEDARKER EAAVAAAAA VPSEPGDPLE AVAFEEKEGK AVLNLLFSR ATKPSALSRA VKVFETFEAK IHHLETRPAQ RPRAGGPHLE YFVRLEVRRG DLAALLSGVR QVSEDEVSPA GPKVPWFPRK VSELDKCHHL VTKFDPDL DL DHPGFSDQVY RQRRKLIAEI AFQYRHGDPI PRVEYTAEEI ATWKEVYTTL KGLYATHACG EHLEAFALLE RFSGYREDNI PQLEDVSRFL KERTGFQLRP VAGLLSARDF LASLAFRVFQ CTQYIRHASS PMHSPEPDCC HELLGHVPML ADRTFAQFSQ DIGLASLGAS DEEIEKLSTL YWFTVEFGLC KQNGEVKAYG AGLLSSYGEL LHCLSEEPEI RAFDPEAAAV QPYQDQTYQS VYFVSESFSD AKDKLRSYAS RIQRPFVSVKF DPYTLAIDVL DSPQAVRRSL EGVQDELDTL AHALSAIG
Source	Mammalian cell
Target Names	TH
Protein Names	Recommended name: Tyrosine 3-monoxygenase EC= 1.14.16.2 Alternative name(s): Tyrosine 3-hydroxylase Short name= TH
Expression Region	2-528
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is involved in the conversion of tyrosine to dopamine. It is the rate-limiting enzyme in the synthesis of catecholamines, hence plays a key role in the physiology of adrenergic neurons. Mutations in this gene have been associated with autosomal recessive Segawa syndrome. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.