



# Recombinant Rat Procathepsin L (CtSl), partial

<b>Product Code</b>	CSB-EP006193RA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P07154
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	IPKTVDWREKGCVTPVKNQGGCGSCWAFSASGCLEGQMFLKTGKLISLSEQN LVDCSHDQGNQGCNGGLMDFAFQYIKENGGDSEESYPYEAKDGSCKYRAE YAVANDTGFVDIPQQEKALMKAVATVGPISVAMDASHPSLQFYSSGIYYEPNC SSKDLDHGVLVVGYYEGT
<b>Source</b>	E.coli
<b>Target Names</b>	CtSl
<b>Protein Names</b>	Recommended name: Cathepsin L1 EC= 3.4.22.15 Alternative name(s): Cyclic protein 2 Short name= CP-2 Major excreted protein Short name= MEPCleaved into the following 3 chains: 1. Procathepsin L 2. Cathepsin L1 heav
<b>Expression Region</b>	114-288aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.