



Recombinant Mouse Creatine kinase M-type (Ckm)

Product Code	CSB-BP005459MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P07310
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	PFGNTHNKF KLNYPQEEY PDLSKHNNHM AKVLTPDLYN KLRDKETPSG FTLDDVIQTG VDNPGHPFIM TVGCVAGDEE SYTVFKDLFD PIIQDRHGGY KPTDKHKTDL NHENLKGGDD LDPNYVLSSR VRTGRSIKGY TLPPhCSRGE RRAVEKLSVE ALNSLTGEFK GKYYPLKSMT EQEQQLIDD HFLFDKPVSP LLLASGMARD WPDARGIWHN DNKSFLVWVN EEDHLRVISM EKGGNMKEVF RRFCVGLQKI EEIFKKAGHP FMWNEHLGYV LTCPSNLGTG LRGGVHVKLA NLSKHPKFEE ILTRLRLQKR GTGGVDTAAV GAVFDISNAD RLGSSSEVEQV QLVVDGVKLM VEMEKKLEKG QSIDDMIPAQ K
Source	Baculovirus
Target Names	Ckm
Protein Names	Recommended name: Creatine kinase M-type EC= 2.7.3.2 Alternative name(s): Creatine kinase M chain M-CK
Expression Region	2-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a cytoplasmic enzyme involved in energy homeostasis and is an important serum marker for myocardial infarction. The encoded protein reversibly catalyzes the transfer of phosphate between ATP and various phosphogens such as creatine phosphate. It acts as a homodimer in striated muscle as well as in other tissues, and as a heterodimer with a similar brain isozyme in heart. The encoded protein is a member of the ATP:guanido phosphotransferase protein family.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.