




Recombinant Dactylis glomeRata Major pollen allergen Dac g 4, partial

Product Code	CSB-YP307868DAC
Abbreviation	Recombinant Dactylis glomerata Major pollen allergen Dac g 4 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P82946
Product Type	Recombinant Proteins
Immunogen Species	Dactylis glomerata (Orchard grass) (Cock's-foot grass)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	DIYNYMEPYVSKVDPTDYFGNEQARTAWVDSGAQLGELSYGVLFNIQYVNYW FAP
Research Area	Others
Source	Yeast
Target Names	N/A
Protein Names	Recommended name: Major pollen allergen Dac g 4 Alternative name(s): Allergen= Dac g 4
Expression Region	1-55aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	8.4kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.