



Recombinant Bacillus sp. Thermostable monoacylglycerol lipase

Product Code	CSB-BP308407BBZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P82597
Product Type	Recombinant Protein
Immunogen Species	Bacillus sp. (strain H-257)
Purity	>85% (SDS-PAGE)
Sequence	SEQYPVLSG AEPFYAENGP VGVLLVHGFT GTPHSMRPLA EAYAKAGYTV CLPRLKGHGT HYEDMERTTF HDWVASVEEG YGWLKQRCQT IFVTGLSMGG TLTLYLAEHH PDICGIVPIN AAVDIPAIAA GMTGGGELPR YLDISGSDLK NPDVKELAYE KTPTASLLQL ARLMAQTAK LDRIVCPALI FVSDVDHVVP PGNADIIFQG ISSTEKEIVR LRNSYHVATL DYDQPMIIE SLEFFAKHAG
Source	Baculovirus
Protein Names	Recommended name: Thermostable monoacylglycerol lipase Short name= MGLP EC= 3.1.1.23
Expression Region	2-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.