



Recombinant Rat Clathrin light chain A (Clta)

| | |
|--------------------------|--|
| Product Code | CSB-EP005591RA |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P08081 |
| Product Type | Recombinant Protein |
| Immunogen Species | Rattus norvegicus (Rat) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MAELDPFGAP AGAPGGPALG NGVAGAGEED PAAFLAQQE SEIAGIENDE AFAILDGGAP GPQAHGEPGP GPDAVDGVMN GEYYQESNGP TDSYAAISEV DRLQSEPESEI RKWREEQTER LEALDANSRK QEAWEKEKAV KELEEWYARQ DEQLQKTKAS NRVADEAFYK QPFADVIGYV TNINHPCYSL EQAAEEAFVN DIDESSPGTE WERVARLCDF NPKSSKQAKD VSRMRSVLIS LKQAPLVH |
| Source | E.coli |
| Target Names | Clta |
| Protein Names | Recommended name: Clathrin light chain A Short name= Lca |
| Expression Region | 1-248 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | Clathrin is a large, soluble protein composed of heavy and light chains. It functions as the main structural component of the lattice-type cytoplasmic face of coated pits and vesicles which entrap specific macromolecules during receptor-mediated endocytosis. This gene encodes one of two clathrin light chain proteins which are believed to function as regulatory elements. Alternative splicing results in multiple transcript variants. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |