



Recombinant Dimethyladenosine transferase 1, mitochondrial (T03F1.7)

Product Code	CSB-MP308594CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P91424
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	IHMY RLRAKKILSQ NYLMDMNITR KIAKHAKVIE KDWVIEIGPG PGGITRAILE AGASRLDVVE IDNRFIPPLQ HLAEAADSRM FIIHQDALRT EIGDIWKNET ARPESVDWHD SNLPAMHVIG NLPFNIASPL IIKYLRDMSY RRGVWQYGRV PLTLTFQLEV AKRLCSPAC DTRSRISIMS QYVAEPKMFV QISGSCFVPR PQVDVGVVRF VPRKTPLVNT SFEVLEKVCR QVFHYRQKYV TKGLKTYLPE ELEDELSDDL LKKCRIDPTT TSIRLGIEQF ADLAEGYNEQ CIRYPGLFLY DYTNKLHNLE DLSKEPNALP PPVPIFAPAP TIDSADNTWS LKNFNCS
Source	Mammalian cell
Target Names	tfbm-1
Protein Names	Recommended name: Dimethyladenosine transferase 1, mitochondrial EC=2.1.1.- Alternative name(s): Mitochondrial 12S rRNA dimethylase 1 Mitochondrial transcription factor B1 S-adenosylmethionine-6-N', N'-adenosyl(rRNA) dimethyltransf
Expression Region	17-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.