



Recombinant Arabidopsis thaliana

Uncharacterized mitochondrial protein AtMg00920 (AtMg00920)

Product Code	CSB-MP308606DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92529
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MKMKSPLEFRG PLVNSSTENP IHILTREKLR HQTSGTNSTE KVILPNILHS YIKNLNLDFT SPYGSTGNNE VLSILRSRFB QSIFPSSGLK CLDTTGDFLI KNVLHKRYES VQKNISNALS SSINSRTAVF FCILFSITVL MEIAPGPLLK KPSLLFSDNL PNVLQYTRDV YVNHVCIHK SLSPCECEEP LNRIIRDMPF QTTFDPLELQ KPSPQ
Source	Mammalian cell
Target Names	AtMg00920
Protein Names	Recommended name: Uncharacterized mitochondrial protein AtMg00920 Alternative name(s): ORF215b
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.