



Recombinant Arabidopsis thaliana Arginine/serine-rich-splicing factor RSP31 (RSP31)

Product Code	CSB-EP308613DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92964
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRPVFVGNFE YETRQSDLER LFDKYGRVDR VDMKSGYAFV YFEDERDAED AIRKLDNFPF GYEKRRLSVE WAKGERGRPR GDAKAPSNLK PTKTLFVINP DPIRTKEHDI EKHFEPYGVK TNVRIRRNFS FVQFETQEDA TKALEATQRS KILDRVVSVE YALKDDDERD DRNGGRSPRR SLSPVYRRRP SPDYGRRPSP GQGRPSPDY GRARSPEYDR YKGPAAYERR RSPDYGRRSS DYGRQRSPGY DRYRSRSPVP RGRP
Source	E.coli
Target Names	RS31
Protein Names	Recommended name: Arginine/serine-rich-splicing factor RSP31
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.