



# Recombinant Bacillus subtilis Peptide deformylase 1 (defA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP308664BRJ   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P94462  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Bacillus subtilis (strain 168)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MAVKKVVTHP AEVLETPAET VTVFDKCLKK LLDDMYDTML EMDGVGLAAP<br>QIGILKRAAV VEIGDDRGRV DLVNPEILEK SGEQTGIEGC LSFNRYGDV<br>TRADYVKVRA FNRQGKPFIL EARGFLARAV QHEMDHLDGV LFTSKISKYY<br>TEDELADMEG   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | defA  |
| <b>Protein Names</b>     | Recommended name: Peptide deformylase 1 Short name= PDF 1 EC= 3.5.1.88<br>Alternative name(s): Polypeptide deformylase 1  |
| <b>Expression Region</b> | 1-160   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |