



Recombinant Methanococcus vannielii Flagellin B2 (flaB2)

Product Code	CSB-EP308692MNS-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P95317
Product Type	Recombinant Protein
Immunogen Species	Methanococcus vannielii (strain ATCC 35089 / DSM 1224 / JCM 13029 / OCM 148 / SB)
Purity	>85% (SDS-PAGE)
Sequence	ASGIGTLI VFIAMVLVAA VAASVLINTS GFLQQKAATT GKESTEQVAS GLQVLQIMGV HDSSNITKLA VYIAPNAGSS AVDLSQSVVT LSNGDVKSIF KYDTTDAAD KSYGAVKTGG SLFAEATLNI SNMTNTQFGI IVIQDADGSC KAATPVLNKG DIVAIVLNLN ETMTEPRTAI TGSIQPEFGA PGIISFTTPA TYLDGSTVVQ LQ
Source	E.coli
Target Names	flaB2
Protein Names	Recommended name: Flagellin B2
Expression Region	13-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.