



Recombinant Schizosaccharomyces pombe Probable 20S rRNA accumulation protein 4 (SPBC25H2.15)

Product Code	CSB-EP308942SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P87156
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSQKVVQTTQA WLGFPDVPIS QEDKPDEYST FLGGFPVFFD GCSLNPNIK CGNCKNLCRL LLQCYAPLEG DNLKERALYV WGCHNPSCR VPNSIVCVRG VRLPLKSDIE AVKSPKAISH LEEKKSSPKE KKVNPFAITS ESSRGLNPFS DATSANPFPS LSTDVNPSKP SSNVFSKPSF AAKAQQSITD QKKTQAKTKT KHIALTTSGM SVHPPVTETY TYPVTEAFQG MFLGLDLEYV PQNKVNSKKD DFSTFKNYTP YLNDSSSEAW KESYEKSPSV YEKTFRLFSE KISHNPTQCL RYERGGTPLL ASGRDKLGQQ LKSVTNFGKS PVPLCPLCKS PRLFEMQLMP HAISILNDEI AEWSTILVAT CSMDCNPPIN KDRVGYAVEW VGIQWD
Source	E.coli
Target Names	SPBC25H2.15
Protein Names	Recommended name: Probable 20S rRNA accumulation protein 4
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.