



Recombinant Cuticle collagen dpy-5 (dpy-5)

Product Code	CSB-BP308977CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P91285
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MVKAVVGFGA ACGISAIVAC LWAALVITND INDMYDDVMG ELGGFRDISD DTWGTLLDVR HGAGESAEQY VRGIFGRHKR SNSQCSCGLP SQGCPAGAPG NPGAPGEPGG TGPDGKNGPT GLPGLNIPIP NDFPKECIKC PAGPPGQDGL PGQEGFQGLP GDAGKRGTPG KDGEPRVGD IGDQGTPGQD GQPGLAGPPG RDGLTGKGQP GVAGRPGMPG PRGEPGNNGN PGEEGQTGAQ GPTGQPGKDG FNGNDGTPGQ AGPQGAVGAD AEYCPCPERK RRRV
Source	Baculovirus
Target Names	dpy-5
Protein Names	Recommended name: Cuticle collagen dpy-5 Alternative name(s): Protein dumpy-5
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.