



# Recombinant Nuclear hormone receptor family member nhr-69 (nhr-69)

<b>Product Code</b>	CSB-YP308983CXY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P91829
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVEEICHICN DKSTGKHYGA ISCDGCKGFF RRSIRKRYHY QCRFEQNCDV TKNKRNACRA CRLQKCVKAG MKSNAIQNER DAIGKRKTT GAEKEDLIDQ LVAAETLCQQ LRSSVIKNTS SLAPYDCGKV KWNYEDARAA TLDDIGKSIH QQLVLFIEWA KSLPQFSFLA QADQAALLKG GAASIIVLGV AYRSICLTVE NTICLANDTL LPKEHATQVG DINCVVGR II DEIVNPMRRL NMDLIEYVAL KAILFFNPVV REINDQSPVE NARYAFLRSL QRRCTDKALE NMEDESMDCR SGKLLLLLPS LQAIAQQLVE DVQLARLFGL VNVDLSMEEL ILNDMKPSDP QILQTSLASP VNSSVKA EVE LEE
<b>Source</b>	Yeast
<b>Target Names</b>	nhr-69
<b>Protein Names</b>	Recommended name: Nuclear hormone receptor family member nhr-69
<b>Expression Region</b>	1-373
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.