



Recombinant Arabidopsis thaliana Arginine/serine-rich-splicing factor RSP41 (RSP41)

Product Code	CSB-MP309000DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92966
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKPVFCGNFE YDARESDLER LFRKYGKVER VDMKAGFAFV YMEDERDAED AIRALDRFEY GRTGRRLRVE WTKNDRGGAG RSGGSRRSSS GLRPSKTLFV INFDAQNTRT RDLERHFEPY GKIVNVRIRR NFAFIQYEAQ EDATRALDAT NSSKLMKVI SVEYAVKDDD SRNGYSPER RRDRSPDRRR RSPSPYRRER GSPDYGRGAS PVAHKRERTS PDYGRGRRSP SPYKRARLSP DYKRDDRRRE RVASPENGAV RNRSPRKGRG ESRSPPPYEK RRESRSPPPY EKRRESRSP PYEKRRERSR SRKSSPENQ QVESPGQIME VEAGRGYDGA DSPIRESPSR SPPAEE
Source	Mammalian cell
Target Names	RS41
Protein Names	Recommended name: Arginine/serine-rich-splicing factor RSP41
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.