



Recombinant Arabidopsis thaliana Putative actin-9 (ACT9)

Product Code	CSB-MP309022DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93738
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKPIVCDKGGH GMVQAGFAGD EAPKVVFPVCV VGRPRDGLNP NESYVGEEGH ANRDILTLKD PMEHGIVNNW DDMEKIWHYT FYNELRVDPE EHPVLLTEAP YNPKANREKM TQIMFESFDV PAMYVSMQSV LYLYSSGRTT GVVLDLGERV SHTVPVYEGY ALPHGILRLD LGGRDLTDYL IEIMTERGYT YT TSAEREIV RDIKEKLCYV ALDYEQEMEK TTKGWTIDKT YVLPDGGQEI T IEAERFMCPE VLFQPSVIGK ESSGIHEATR NSILKCPVDT RRD MYGNILM TGGTTMLHGI KERMTKELNA LVPSSMKVKV VVPPSECSV WIGGSILASL STFHQMWITK DEYEEHGA AI VHRKCV
Source	Mammalian cell
Target Names	ACT9
Protein Names	Recommended name: Putative actin-9
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.