



Recombinant Arabidopsis thaliana 1-aminocyclopropane-1-carboxylate oxidase homolog 6 (At1g04350)

Product Code	CSB-MP309024DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93824
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	METKEFDSYS ERKAFDETKT GVKGLIDAH I TEIPRIFCLP QGSLSDKKPF VSTTDFAIPI IDFEGLHVSR EDIVGKIKDA ASNWGFQVI NHGVPLNVLQ EIQDGVRRFH EEAPEVKKTY FTRDATKRFV YNSNFDLYSS SSCVNWRDSF ACYMAPDPPN PEDLPVACRV AMFEYSKHMM RLGDLLFELL SEALGLRSDK LKSMDCMKGL LLLCHYYPPC PQLDLTIGTN NHSDNSFLTI LLQDQIGGLQ IFHQDCWVDV SPIPGALVIN MGDFLQLITN DKVISVEHRV LANRAATPRI SVASFFSTSM RPNSTVYGPI KELLSEENPS KYRVIDLKEY TEGYFKKGLD GTSYLSHYKI
Source	Mammalian cell
Target Names	At1g04350
Protein Names	Recommended name: 1-aminocyclopropane-1-carboxylate oxidase homolog 6 EC= 1.14.-.-
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.