



Recombinant Human Nuclear factor 1 C-type (NFIC)

Product Code	CSB-BP015756HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P08651
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MYSSPLCLTQ DEFHPFIEAL LPHVRAFAYT WFNLQARKRK YFKKHEKRMS KDEERAVKDE LLGEKPEVKQ KWASRLAKL RKDIRPECRE DFVLSITGKK APGCVLSNPD QKGKMRRIDC LRQADKVVRL DLVMVILFKG IPLESTDGER LVKAAQCGHP VLCVQPHHIG VAVKELDLYL AYFVRERDAE QSGSPRTGMG SDQEDSKPIT LDTTDFQESF VTSGVFSVTE LIQVSRTPVV TGTGPNFSLG ELQGHLAYDL NPASTGLRRT LPSTSSSGSK RHKSGSMEED VDTSPGGDYY TSPSSPTSSS RNWTEDEMEGG ISSPVKKTEM DKSPFNPSPP QDSPRLSSFT QHHRPVIAVH SGIARSPHPS SALHFPTTSI LPQTASTYFP HTAIRYPHPL NPQDPLKDLV SLACDPASQQ PGPLNGSGQL KMPSHCLSAQ MLAPPPGLP RLALPPATKP ATTSEGGATS PTSPSYSPPD TSPANRSFVG LGPRDPAGIY QAQSWYLG </p>
Source	Baculovirus
Target Names	NFIC
Protein Names	Recommended name: Nuclear factor 1 C-type Short name= NF1-C Short name= Nuclear factor 1/C Alternative name(s): CCAAT-box-binding transcription factor Short name= CTF Nuclear factor I/C Short name= NF-I/C Short
Expression Region	1-508
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.