



Recombinant Treponema pallidum Periplasmic zinc-binding protein troA (troA)

Product Code	CSB-MP309098TPR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P96116
Product Type	Recombinant Protein
Immunogen Species	Treponema pallidum (strain Nichols)
Purity	>85% (SDS-PAGE)
Sequence	FGSKDAAA DGKPLVVTI GMIADAVKNI AQGDVHLKGL MGPGVDPHLY TATAGDVEWL GNADLILYNG LHLETKMGEV FSKLRGSRLV VAVSETIPVS QRLSLEEAEF DPHVWFDVKL WSYSVKAVYE SLCKLLPGKT REFTQRYQAY QQQLDKLDAY VRRKAQSLPA ERRVLVTAHD AFGYFSRAYG FEVKGLQGVS TASEASAHDM QELAAFIAQR KLPAIFIESS IPHKNVEALR DAVQARGHVV QIGGELFSDA MGDAGTSEGT YVGMVTHNID TIVAALAR
Source	Mammalian cell
Target Names	troA
Protein Names	Recommended name: Periplasmic zinc-binding protein troA Alternative name(s): Tromp-1
Expression Region	23-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.