



Recombinant Rabbit Cytochrome P450 2E1 (CYP2E1)

Product Code	CSB-MP006425RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P08682
Product Type	Recombinant Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	≥85% (SDS-PAGE)
Sequence	AVLGITVAL LGWMVILLFI SVWKQIHSSW NLPPGPFPLP IIGNLLQLDL KDIPKSFGR L AERFGPVFTV YLGSRRVVVL HGYKAVREML LNHKNEFSGR GEIPAFREFK DKGIIFNNGP TWKDRRFS L TTLRDYGMGK QGNEDRIQKE AHFLLLELRK TQGQPFDP TF VIGCTPFNVI AKILFNDRFD YKDKQALRLM SLFNENFYLL STPWLQVYNN FSNYLQYMPG SHRKVIKNVS EIKEYTLARV KEHHKSLDPS CPRDFIDSL L IEMEKDKHST EPLYTLENIA VTVADMFFAG TETTSTTLRY GLLILLKHPE IEEKLHEEID RVIGPSRMP S VRDRVQMPYM DAVVHEIQRF IDLVPSNLPH EATRDTTFQG YVIPKGT VVI PTLDSL LLYDK QEFDPPEKFK PEHFLNEEGK FKYS DYFKPF SAGKRVCVGE GLARME LFL L LSAILQHFN L KPLVDPEDID LRNITVGFGR VPPRYKLCVI PRS
Source	Mammalian cell
Target Names	CYP2E1
Protein Names	Recommended name: Cytochrome P450 2E1 EC= 1.14.13.- Alternative name(s): 4-nitrophenol 2-hydroxylase EC= 1.14.13.n7 CYP11E1 Cytochrome P450 isozyme 3A Short name= Cytochrome P450 LM3A Cytochrome P450-ALC
Expression Region	2-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and is induced by ethanol, the diabetic state, and starvation. The enzyme metabolizes both endogenous substrates, such as ethanol, acetone, and acetal, as well as exogenous substrates including benzene, carbon tetrachloride, ethylene glycol, and nitrosamines which are premutagens found in cigarette smoke. Due to its many substrates, this enzyme may be involved in such varied processes as gluconeogenesis, hepatic cirrhosis, diabetes, and cancer.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.