



# Recombinant Tubulin alpha-3 chain (mec-12)

<b>Product Code</b>	CSB-MP309355CXY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P91910
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MREVISIHIG QAGVQIGNAC WELYCLEHGI QPDGQMPSDK SLGGSDDSFSTFFSETGSGR HVPRAVMVDL EPTVIDEIRT GTYRSLFHPE QLITGKEDAA NNYARGHYTI GKEIDLTLDRIRRLADNCT GLQGFVLFHS FGGGTGSGFT SLLMERLSVD YGKKAKLEFS IYPAPQVSTA VVEPYNSILT THITTLEHSDC SFMVDNEAIY DICRRNLIDIE RPSYTNLNRL IGQIVSSITA SLRFDGALNV DLTEFQTNLV PYPRIHFPLA TFSPVISA EK AYHEQLSVAE ITNMCFEPHN QMVKCDPRHG KYMAVCLLFR GDVVPKDVNA AIATIKTKRS IQFVDWCPTG FKVGINYQPP TVVPGGDLAK VPRAVCMLSN TTAIAEAWAR LDHKFDLMYA KRAVHWYVG EGMEEGEFSE AREDLALEK DYEEVGVD SM EDNGEEGDEY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	mec-12
<b>Protein Names</b>	Recommended name: Tubulin alpha-3 chain Alternative name(s): Mechanosensory abnormality protein 12 Short name= MEC-12
<b>Expression Region</b>	1-450
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.