



# Recombinant Human Homeobox protein Hox-C4 (HOXC4)

<b>Product Code</b>	CSB-MP010675HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P09017
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIMSSYLMDS NYIDPKFPPC EEYSQNSYIP EHSPEYYGRT RESGFQHHHQ ELYPPPPRP SYPERQYSCT SLQGPGNSRG HGPAQAGHHH PEKSQSLCEP APLSGASASP SPAPPACSQP APDHPSSAAS KQPIVYPWMK KIHVSTVNPN YNGGEPKRSR TAYTRQQVLE LEKEFHYNRY LTRRRRIEIA HSLCLSERQI KIWQNRMRK WKKDHRLPNT KVRSAPPAGA APSTLSAATP GTSEDHSQSA TPPEQQRAED ITRL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HOXC4
<b>Protein Names</b>	Recommended name: Homeobox protein Hox-C4 Alternative name(s): Homeobox protein CP19 Homeobox protein Hox-3E
<b>Expression Region</b>	1-264
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene, HOXC4, is one of several homeobox HOXC genes located in a cluster on chromosome 12. Three genes, HOXC5, HOXC4 and HOXC6, share a 5 non-coding exon. Transcripts may include the shared exon spliced to the gene-specific exons, or they may include only the gene-specific exons. Two alternatively spliced variants that encode the same protein have been described for HOXC4. Transcript variant one includes the shared exon, and transcript variant two includes only gene-specific exons.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.