



Recombinant Arabidopsis thaliana

Uncharacterized mitochondrial protein AtMg00810 (AtMg00810)

Product Code	CSB-BP309363DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92519
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MYLLLYVDDI LLTGSSNTLL NMLIFQLSST FSMKDLGPVH YFLGIQIKTH PSGLFLSQTK YAEQILNAG MLDCKPMSTP LPLKLNSSVS TAKYPDPSDF RSIVGALQYL TLTRPDISYA VNIVCQRMHE PTLADFDLLK RVLRYVKGTI FHGLYIHKNS KLNQAFCDSDWAGCTSTRR STTGFCFTFLG CNIISWSAKR QPTVSRSTE TEYRALALTA AELTWSSASR SRDPSAMNTN
Source	Baculovirus
Target Names	AtMg00810;
Protein Names	Recommended name: Uncharacterized mitochondrial protein AtMg00810 Alternative name(s): ORF240b
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.