



Recombinant Arabidopsis thaliana Germin-like protein subfamily 3 member 1 (GLP1)

Product Code	CSB-BP309401DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P94040
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	VQ DFCVANLKRA ETPAGYPCIR PIHVKATDFV FSGLGTPGNT TNIINAAVTP AFAAQFPGLN GLGLSTARLD LAPKGVIPMH THPGASEVLF VLTGSITAGF VSSANAVYVQ TLKPGQVMVF PQGLLHFQIN AGKSSASAVV TFNSANPGLQ ILDFALFANS LPTELVVGTT FLDATTVKKL KGVLGGTG
Source	Baculovirus
Target Names	GLP1
Protein Names	Recommended name: Germin-like protein subfamily 3 member 1 Short name= At-GERM1 Short name= AtGER1 Short name= AtGLP1
Expression Region	19-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.