



# Recombinant Mouse Cholecystokinin (Cck)

<b>Product Code</b>	CSB-EP004771MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P09240
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	QPVVPAEATD PVEQRAQEAP RRQLRAVLRT DGEPRARLGA LLARYIQQVR KAPSGRMSVL KNLQSLDPSH RISDRDYMGW MDFGRRSAED YEYPS
<b>Source</b>	E.coli
<b>Target Names</b>	Cck
<b>Protein Names</b>	Recommended name: Cholecystokinin Short name= CCK Cleaved into the following 3 chains: 1. Cholecystokinin-33 Short name= 2. CCK33 3. Cholecystokinin-12 Short name= 4. CCK12 5. Cholecystokinin-8 Short name= 6. CC
<b>Expression Region</b>	21-115
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Cholecystokinin is a brain/gut peptide. In the gut, it induces the release of pancreatic enzymes and the contraction of the gallbladder. In the brain, its physiologic role is unclear. The cholecystokinin pro-hormone is processed by endo- and exo-proteolytic cleavages.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.