



Recombinant Mouse Phosphoglycerate kinase 1 (Pgk1)

Product Code	CSB-BP017856MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P09411
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SLSNKLTL D KLDVKGKRVV MRVDFNVPMK NNQITNNQRI KAAVPSIKFC LDNGAKSVVL MSHLGRPDGV PMPDKYSLEP VAAELKSLLG KDVLFLKDCV GPEVENACAN PAAGTVILLE NLRFHVEEEG KGKDASGNKV KAEPKIDAF RASLSKLGDV YVNDAFGTAH RAHSSMVGVN LPQKAGGFLM K KELNYFAKA LESP RPFLA ILGGAKVADK IQLINNMLDK VNEMIIGGGM AFTFLKVLNN MEIGTSLYDE EGAKIVKDLM SKAEKNGVKI TLPVDFVTAD KF DENAKTGQ ATVASGIPAG WMGLDCGTES SKKYAEAVGR AKQIVWNGPV GVFEWEAFAR GTKSLMDEVV KATSRGCITI IGGGDTATCC AKWNTEDKVS HVSTGGGASL ELLEGKVLPG VDALS NV
Source	Baculovirus
Target Names	Pgk1
Protein Names	Recommended name: Phosphoglycerate kinase 1 EC= 2.7.2.3
Expression Region	2-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a glycolytic enzyme that catalyzes the conversion of 1,3-diphosphoglycerate to 3-phosphoglycerate. The encoded protein may also act as a cofactor for polymerase alpha. This gene lies on the X-chromosome, while a related pseudogene also has been found on the X-chromosome and another on chromosome 19.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.