



Recombinant 39S ribosomal protein L41, mitochondrial (mrpl-41)

Product Code	CSB-YP309716CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P90993
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	PWPFVKKGA FGQRKIGPMN YEKFKWPGQN REFPELSPKF QKLNPKELHR YTGTVQADGFH DEKTGEFVPV KEMRSELVVP NLDGFKLRPY VSYRTDVQIE KRRVAYEKKV LEKGSERLAD LHTVEDERWP PPKMSAETLF ELAYGDTVRS AYKEGKYGNK
Source	Yeast
Target Names	mrpl-41
Protein Names	Recommended name: 39S ribosomal protein L41, mitochondrial Short name= L41mt Short name= MRP-L41
Expression Region	22-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.