



Recombinant Bacillus subtilis Uncharacterized protein yeaD (yeaD)

Product Code	CSB-YP309792BRJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P94475
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	>85% (SDS-PAGE)
Sequence	MNKKWQWLVIY GWKLIILTLL TAAVFSYAMF QGGFVSWFLF YAFLPFVLYA GLLALYPLRS FQASRQMDKT QLHAGDRLGV TVTLRRKLPF PLMYMVIEDC LPEGIESLNR DGAAAKRLVF PWFKRSMTMS YELARVPRGE HHFHSVRVRT GDVLGLIEKT AFFELDDTLF VYPFYQRFSY QVNERHQEDG VSGSSPIHQH HSSVAASVRN YQPGDRFAAL DWKTSARRSQ LMTKEFEPSR SKNLFLLMDR FSSDAFEEVV SVTASILHSV LKNGAGAGLA SIGKEKNIFP IQEGDQHFKH MLRHIAIAHC DAADPISRYA REELGKPSVR QADQVVVTGQ LTEDMLHLAE IGGGRVTVIL AKEKDAELSQ AENVMIERMM KRQIRVRIMR GGRVSRVV
Source	Yeast
Target Names	yeaD
Protein Names	Recommended name: Uncharacterized protein yeaD
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.