



Recombinant Rat Glutamine synthetase (Glul)

Product Code	CSB-EP009553RA-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P09606
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	ATSASSHLN KGIKQMYMNL PQGEKIQLMY IWVDGTGEGEGL RCKTRTLDCD PKCVEELPEW NFDGSSTFQS EGSNSDMYLH PVAMFRDPFR RDPNKLVFCE VFKYNRKPAE TNLRHSCRKRI MDMVSSQHPW FGMEQEYTLM GTDGHPFGWP SNGFPGPQGP YYCGVGADKA YGRDIVEAHY RACLYAGIKI TGTNAEVMPA QWEFQIGPCE GIRMGDHLWV ARFILHRVCE DFGVIATFDP KPIPGNWNGA GCHTNFSTKA MREENGLRCI EEAIKLSKR HQYHIRAYDP KGGLDNARRL TGFHETSNIN DFSAGVANRS ASIRIPRIVG QEKKGYFEDR RPSANCDPYA VTEAIVRTCL LNETGDEPFQ YKN
Source	E.coli
Target Names	Glul
Protein Names	Recommended name: Glutamine synthetase Short name= GS EC= 6.3.1.2 Alternative name(s): Glutamate decarboxylase EC= 4.1.1.15 Glutamate-- ammonia ligase
Expression Region	2-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.