



# Recombinant Chicken Vimentin (VIM)

<b>Product Code</b>	CSB-EP025857CH-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P09654
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SFTSSKNSS YRRMFGGGSR PSSGTRYITS STRYSLGSAL RPSSARYVSA SPGGVYRTKA TSVRLRSSMP PMRMHDAVDF TLADAINTEF KANRTNEKVE LQELNDRFAN YIDKVRFLAQ QNKILLAELE QLKGGKTSRL GDLYEEEMRD VRRQVDQLTN DKARVEVERD NLADDIMRLR EKLQEEMLQR EEAESTLQSF RQDVDNASLA GLDLERPVES LQEEIVFLKK LHDEEIRELQ AQLQEHIQI DMDVSKPDLT AALRDVRQQY ESVAAKNLQE AEEWYKSKFA DLSEAANRNN DALRQAKQEA NEYRRQIQSL TCEVDALKGS NESLERQMRE MEENFAVEAA NYQDTIGRLQ DEIQNMKEEM ARHLREYQDL LNVKMALDIE IATYRKLLG EESRINMPIP TFASLNLRET NIESQPIVDT HSKRTLLIKT VETRDGQVIN ETSQHHDDLE
<b>Source</b>	E.coli
<b>Target Names</b>	VIM
<b>Protein Names</b>	Recommended name: Vimentin
<b>Expression Region</b>	2-460
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.