



Recombinant Danio rerio Borealin-2 (cdca9)

Product Code	CSB-MP309976DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P86347
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAPRRTRKVS QDSDGQADDQ HSFEQKIRLT KRKELFIQQF EKEAQDRINE MEANLNKLLA TVDRVFKIEL MKMPLSLHTT LIKDLMNDDD TSVG EVT MAL KCASPEIQKP LSRKPSKKAL NALAGQQRSS SQSKTPIEGQ KKPTKKT LHS SKSTGSLRCA STINAKRTQG RVVKLSDQAN ALGVQFRQTS RSVGDEL MMA TATIVTSHGE TLFLSEDNKD EINVELLDDA AVNQMRKIKE LMDYLCNKVR INNTC
Source	Mammalian cell
Target Names	cdca9
Protein Names	Recommended name: Borealin-2 Alternative name(s): Cell division cycle-associated protein 8.2 Cell division cycle-associated protein 9 Dasra-A Short name= DrDasraA
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.