



Recombinant Schizosaccharomyces pombe Pre-mRNA-splicing factor cwc2 (cwf2)

Product Code	CSB-MP310043SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P87126
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSENGLEQEV TVEEKNNVDVT EKILVEGEKS KEYEETPRKV KIVKRKKQPA RKQIETRPEY EMEPEQPGQV YNLWYNKWSG GMRQDPLKSQ VKSETRCVIS RDSGYTKADK NPGSFFCLYF ARGMCSEGSK CEYLHRLPKD TDDFNANVDC FGREKHADYR DDMGGVGSFL RQNYTLYVGG IPTDDIEEI VSRHFAEWGD IERIRVLNSR GIAFITYLNE ANAQFAKEAM AHQSLDHDEC LNVRWATTDTP NPASQARNQR RLEERAANAV KKLLPKQFLL DLEETKNGKS GNRKRKLELE FGLKGYVPSD DLLYADGANS VHNQLAANEF PNKSQSEEGS NDDHKSVT ESQNKVNSQ ILSDLQVAKQ AVHTNQSALV SYYDSDED
Source	Mammalian cell
Target Names	cwf2
Protein Names	Recommended name: Pre-mRNA-splicing factor cwc2 Alternative name(s): Cell cycle control protein cwf2 Complexed with cdc5 protein 2 Pre-mRNA-processing protein 3
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.