



Recombinant Putative thiosulfate sulfurtransferase F11G11.9 (F11G11.9)

Product Code	CSB-MP310080CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P91247
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MAVDMIVEPK WVVQNFNGNIR ILDAWTFKP KADVAEYKAK YYNKFGVGMN ELKNPEYLAE HINGAAHFNF DIAYYPSEDE RFTLYTPEEF SSVVKRLGVF NGDHLVIYGR GKDGGMAAAS RAYWTFRYYG YTTVSVLNGG IEAFKLAQGV VQSDSKAEGI RCKDAIHFPI GEVCAAKGFK KKTDCDQAF A AKGIKVGDTV VISCIGLSA SAICLAAARS GIVAKLYNGG VHELAYKAPQ HLNMRVTL L L HAITVLRCTY IHFTFLYAIV IKIERIV
Source	Mammalian cell
Target Names	mpst-4
Protein Names	Recommended name: Putative thiosulfate sulfurtransferase F11G11.9 EC= 2.8.1.1
Expression Region	1-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.