



Recombinant Probable phosphoserine aminotransferase (F26H9.5)

Product Code	CSB-EP310087CXY-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P91856
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MAAPGRKINF AAGPAKLPEE VLLKMQEEQL NFNNLGVSVI EMSHRSKEFG ALLNETISLI RELMNVPDNF EILFMQGGGT GQFAAIPLNL KGDHEHADI VTGAWSSKAA DEAGKYINVK KVFQPSKPYV TVPDQENWVH DEKAAAYLYYC ANETVHGIEF TPTAPESHNV PLVADVSSNF MARPFDKDH GVVFGGAQKN LGAAGLTIVI VRKDLIGKQQ AITPSVFSYK EMIANNSLYN TPPTGGIYTT NLVLKWIKSK GGLQAIYELN LQKSGMIYDI IDNSNGFYHC AVDKRYRSIM NVCFRIGGPS GNDELEEKFL KGSIERNMIS LKGHRVGGI RASLYNAISV EETQVLATWM NEFQKLHNTN
Source	E.coli
Target Names	F26H9.5
Protein Names	Recommended name: Probable phosphoserine aminotransferase Short name= PSAT EC= 2.6.1.52 Alternative name(s): Phosphohydroxythreonine aminotransferase
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.